



VICTOR - G. & COMPANY

Super - G. Stainless Steel HATCHERIES Needles



COMPONENTS:

HUBS

Made from free cutting brass as per IS:319-1974, is properly electroplated to make it corrosion resistant and suitable for Luer Lock or Luer Tip glass syringe for LEAK PROOF JOINT (as per IS:3234-1965). Hubs with RECORD FITTING ARE ALSO AVAILABLE

Hatcheries Needles feature over- sized Chrome Plated Hubs, for secure gripping even under the most difficult circumstances. Needles with standard (regular) hubs are also available in all sizes configuration, as they are useful when working with smaller animals .

CANNULA

Made from hyperchrome stainless steel tubing, our cannulae have super sharp LANCENT POINT, offering longer life, reduced coring and painless penetration. Points are machine ground to assure uniform sharpness and minimal patient discomfort due to tearing of skin. Precision grinding, micro bead jet blasting and ultrasonic cleaning assure freedom from burrs and debris.

SS tubing is fabricated without exposing the interior to lubricants or mechanical drawing, offering, therefore, the cleanest, purest tube interiors available anywhere.

PACKING

12 needles are packed in Plastic Box. 12 Boxes in export worthy duplex box.



STILLET / CLEANING WIRES are supplied on request

Poultry & Hatcheries (birds) needles are available in 18G (1.3mm) to 23G (0.6mm) x 1/4"(6mm) - 1/2" (13mm) with standard luer lock hub, and standard bevel cannula.

These stainless steel needles are designed to maintain their sharpness through a large number of injections. There are a wide range of sizes available depending upon your specific requirements. They can be relied upon to provide the most consistent sharpness and length of needle.

Vaccination is usually done with a reusable stainless steel needle attached to a multiple injector. The number of small animals that can be injected with each needle is 500 - 1000. This number depends upon the skill of the vaccination team and the size of the animal. Care should be taken to avoid dulling the needle by injecting through scales. Any accumulated scales and debris should be carefully removed from the needle to maximize needle life and avoid reducing the penetration depth. The length of needle chosen will vary depending upon the species being vaccinated. Needle diameter should be as small as possible while allowing the vaccine to pass through easily.

The needles manufacture process totally eliminates interior coarseness of the needle tubing making it an excellent choice for critical blood flow applications and delicate procedures where turbulence and cell damage must be minimized.



47, INDUSTRIAL ESTATE, KANPUR 208 012, INDIA

TEL + 91 512 2221481, 2672009 FAX + 91 512 2217250

Email info@victor-g.com, victorg@sancharnet.in

Website www.victor-g.com, www.veterinaryneedles.com



VICTOR - G. & COMPANY

Super - G. Stainless Steel HATCHERIES Needles

PRODUCT RANGE

6G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
8G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
10G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
11G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
12G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
13G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
14G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
15G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
16G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
17G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
18G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
19G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
20G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
21G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
22G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
23G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
24G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
25G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
26G	1/4"	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"

Veterinary Needles
Luer Lock
Standard Sizes

6G	From 2" upto 12"
8G	From 2" upto 12"
10G	From 2" upto 12"
11G	From 2" upto 12"
12G	From 2" upto 12"
13G	From 2" upto 12"
14G	From 2" upto 12"
15G	From 2" upto 12"
16G	From 2" upto 12"
17G	From 2" upto 12"
18G	From 2" upto 12"
19G	From 2" upto 12"
20G	From 2" upto 12"
21G	From 2" upto 12"
22G	From 2" upto 12"
23G	From 2" upto 12"
24G	From 2" upto 12"
25G	From 2" upto 12"
26G	From 2" upto 12"

Veterinary Needles
Luer Lock
Long Length

POINT STYLES

Long Bevel



Intra-Muscular
12° - 14°

Short Bevel



Intravenous
15° - 19°

Intradermal



Intradermal
30°

Flat Point 45°



Tubing
Connector

Blunt 90°



Syringe
Pipetting

Bent End



with Blunt
Point

Pencil Point




Spinal

Deflected Point



Huber

Plugged Point



Irrigation

SPECIALITY NEEDLES

HYPODERMIC
VETERINARY
SPINAL
CUSTOM FABRICATED
LARGE HUBS
LUER LOCK ADAPTERS
BLUNT END
PENCIL POINT
HYPODERMIC NEEDLE TUBING

Needles of special lengths and sizes can be manufactured on demand.
Call / Fax / Email your special requirements. ALL NEEDLES FURNISHED NON STERILE

47, INDUSTRIAL ESTATE, KANPUR 208 012, INDIA

TEL + 91 512 2221481, 2672009 FAX + 91 512 2217250

Email info@victor-g.com, victorg@sancharnet.in

Website www.victor-g.com, www.veterinaryneedles.com